

PRODUCT INFORMATION

Common Name	D2E7
Synonyms	DIF;TNF-alpha;TNFA;TNFSF2;cachexin;cachectin;TNF α
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	TNFA
Uniprot ID	P01375
Description	PE-conjugated Anti-TNF α (Adalimumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

