

PRODUCT INFORMATION

Common Name	BSI-045B;BSI-04502;BSI045B;BSI04502;TQC2731;TQC-2731
Synonyms	TSLP
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TSLP
Uniprot ID	Q969D9
Description	PE-conjugated Anti-TSLP(bosakitug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

