

**PRODUCT INFORMATION**

<b>Common Name</b>	AMG 157, AMG-157, MEDI-9929, MEDI9929, anti-TSLP-Receptor (TSLP-R; CrI2), tezepelumab-ekko
<b>Synonyms</b>	Thymic stromal lymphopoietin
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Homo sapiens
<b>IgG type</b>	IgG2
<b>Reactivity</b>	Human
<b>Target</b>	TSLP
<b>Uniprot ID</b>	Q969D9
<b>Description</b>	PE-conjugated Anti-TSLP(tezepelumab biosimilar) mAb
<b>Delivery</b>	3-4 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

DIMABIO CONFIDENTIAL

