

PRODUCT INFORMATION

Common Name	SSS-23;SSS23;SSS 23;TTAC-0001;TTAC0001;TTAC 0001;TTAC-1;TTAC1;TTAC 1
Synonyms	FLK1;CD309;VEGFR;VEGFR2
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	VEGFR2
Uniprot ID	P35968
Description	PE-conjugated Anti-VEGFR2(olinvacimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

